

Suspension Hepatocyte Stability Incubations Protocol

Equipment and Reagents needed:

- Standard tissue culture treated 24-well plates
- Micropipette set
- Multi-channel pipettor
- Refrigerated bench top centrifuge
- 37°C Incubator
- Hepatocytes: ≥1 vial based on assay needs
 - At least 5.0×10^6 cells are needed per compound, assuming 5 time points, +controls and viability checks
- Hepatocyte incubation medium, depending on assay needs:
 - INVITROGRO CP (contains serum)
 - INVITROGRO HT (contains serum)
 - INVITROGRO HI (serum free)
 - INVITROGRO KHB (serum free)
 - CyrostaX Hepatocyte Incubation Medium (serum free)
- Test compound(s) and solvent
- Stop solution

Experimental Design and Technical Tips:

Purpose: to monitor test compound stability over time using human primary hepatocytes as a metabolic matrix.

- Typical incubation reactions contain 0.25×10^6 cells per well (250µL total) in a 24-well format.
- Common time points tested include 0, 10, 60, 120, and 240 min. (the 0-time point, non-incubated is 100% of starting compound).
 - Each time point should use a dedicated plate.
 - 3 replicate wells per time point.
- Positive controls: 7-Ethoxycoumarin (1.0 µM) or Midazolam (1.0 µM) in duplicate at 0 and 60 min.
- Negative control is vehicle only.
- Make sure the organic solvent used for test compound stocks is ≤ 0.1% total volume.
- Stop reaction, acetonitrile is commonly used as the stop solution.
 - Spike stop solution with internal standard for Metabolite ID (14C 7-ethoxycoumarin).
- Spin samples at 920 x g for 10 minutes at 10°C to separate protein precipitates from solubilized metabolite(s).
- For analytical assessment, there will be 19 wells/tubes per reaction, plus viability checks at each time point (zero substrate samples).

General Protocol

1. Warm 10mLs of incubation medium to 37°C
 - a. This is enough medium for 1 test compound with 5 time points.
2. Prepare an appropriate volume of incubation media with test compound and keep at RT.

General Protocol Continued

3. Pipette 150 μ L medium containing test compound into wells.
 - a. For other well plate formats, 3/5 of the total well volume will be incubation media and test article.
4. Incubate plate at 37°C for 15 minutes.
5. Add stop reagent to 0 time point wells/plate (250 μ L or 1 volume).
6. Thaw, centrifuge, and dilute hepatocytes to 2.5×10^6 /mL in sterile, warmed incubation media.
 - a. Keep cells in 37°C incubator until needed, cells can be thawed prior to adding test article into plates.
7. Dispense 100 μ L of hepatocytes into each well for a final concentration of 0.25×10^6 hepatocytes/well.
 - a. For other well plate formats, 2/5 of the total well volume will consist of dilute hepatocytes.
8. Harvest 0 time point by adding cells to stop solution and collect samples.
9. Move plate(s) into 37°C culture incubator
 - a. When the medium is bicarbonate buffered for non-atmospheric CO₂ concentrations, a CO₂ incubator (typically 5%) should be used.
10. Rock plates at 150 RPM during incubation.
11. At appropriate times, add 1 volume of stop solution and collect the contents of the wells.
12. Spin down samples at 920 x g for 10 minutes at 10°C and run supernatant on MS/MS.

Viability wells are on a separate plate with only incubation media; no test article is added unless toxicity data is also collected at the same time.