



Product Information Sheet

Rat CYP2C11LR Bactosomes

Catalogue Number: CYP365
Product Description: Rat CYP2C11 and rat CYP-reductase coexpressed in *Escherichia coli*
Amount: 1 nmol cytochrome P450
Storage Buffer: 50 mM Tris-acetate (pH 7.6), 250 mM sucrose, 0.25 mM EDTA
Storage: Store at -80°C. Avoid frequent temperature changes. Thaw on ice.
Intended Use: For laboratory (research) purposes only.

Batch Number: RC2C11LR001
P450 Concentration: 2.5 nmol/ml
Protein Concentration: 12.9 mg/ml
Specific P450 Content: 194 pmol/mg protein
Cytochrome c Reductase Activity: 3054 nmol/min/mg protein

¹P450 Activity Data

Activity: Testosterone 16 α -hydroxylase
²Linear with time up to: 12 min
³Linear with P450 up to: 20 pmol/ml
 V_{max} : 28 pmol/min/pmol CYP
 K_m : 6.6 μ M

Approved by 

Date 22nd June 2022

Notes:

¹Unless otherwise stated, all assays are carried out at 37°C in 100 mM potassium phosphate buffer (pH 7.4) containing 5 mM magnesium chloride and an NADPH generating system [5 mM glucose 6-phosphate, 1 mM NADP⁺, 1 U/ml glucose 6-phosphate dehydrogenase].

²Determined at 5 μ M substrate, using a P450 concentration of 6.25 pmol/ml.

³Determined at 5 μ M substrate, using an incubation time of 5 min.

Safety Data

The toxicological properties of this reagent have not been investigated. Exercise due care when handling.

IMPORTANT DISCLAIMER

1. Products supplied by Cypex may be harmful if misused. Any product ordered from Cypex must not be used for any purpose other than the intended use specified on the Product Information Sheet which accompanies the product when delivered. Please ensure that products supplied by Cypex are used safely, and, in particular, that the products do not come into direct human contact.

2. Cypex undertakes that all products supplied by it shall conform to the specification appearing on the Product Information Sheet which accompanies the product when delivered. In the event that products supplied by Cypex fail to meet such specification, Cypex will supply replacement products free of charge provided that the defective products have been returned to Cypex within a reasonable time of delivery (Cypex will reimburse reasonable costs incurred in returning defective products), and this shall be the only remedy available in respect of the supply of any defective products.

3. Cypex gives no undertakings or warranties in respect of products supplied other than the undertaking given at 2, above. Furthermore, Cypex has given and gives no undertaking or warranty nor does it make any representation, directly or indirectly, express or implied, about the quality, character or suitability of products supplied by Cypex.

4. Cypex shall have no liability in respect of loss or damage (whether direct or indirect or in the nature of loss of profit) suffered as a result of the supply of products by Cypex, or the subsequent use, storage or other application of the products, whether or not the products meet the specification referred to at 2, above.