

oYo-Link® Antibody Labeling Products

Antibody-Drug Conjugation Reagents

Conjugate drugs directly to your antibody of choice. Produce conjugates with defined antibody-drug ratios.

- ✓ High potency ADCs
- ✓ Rapid ADC prototyping: 30 seconds hands on time
- ✓ Site-specific labeling ensures homogenous ADCs



oYo-Link® VcMMAE

oYo-Link VcMMAE is an antibody-drug conjugation reagent enabling site-specific labeling of compatible antibodies with up to two VcMMAE drugs per antibody for subsequent cell killing assays or other RUO drug delivery applications.



oYo-Link® DM1

oYo-Link DM1 is an antibody-drug conjugation reagent enabling site-specific labeling of compatible antibodies with up to two DM1 drugs per antibody for subsequent cell killing assays or other RUO drug delivery applications.



oYo-Link® VcMMAF

oYo-Link VcMMAF is an antibody-drug conjugation reagent enabling site-specific labeling of compatible antibodies with up to two VcMMAF drugs per antibody for subsequent cell killing assays or other RUO drug delivery applications.



oYo-Link® DM4

oYo-Link DM4 is an antibody-drug conjugation reagent enabling site-specific labeling of compatible antibodies with up to two DM4 drugs per antibody for subsequent cell killing assays or other RUO drug delivery applications.

Product	Catalog #	Amount Labeled
VcMMAE	AT7001-100	Label 100 µg Antibody
VcMMAE	AT7001-mIgG1-100	Label 100 µg mouse IgG1 Antibody
VcMMAE	AT7001-500	Label 500 µg Antibody
VcMMAE	AT7001-mIgG1-500	Label 500 µg mouse IgG1 Antibody
VcMMAE	AT7001-1000	Label 1000 µg Antibody
DM1	AT7002-100	Label 100 µg Antibody
DM1	AT7002-mIgG1-100	Label 100 µg mouse IgG1 Antibody
DM1	AT7002-500	Label 500 µg Antibody
DM1	AT7002-mIgG1-500	Label 500 µg mouse IgG1 Antibody
DM1	AT7002-1000	Label 1000 µg Antibody
VcMMAF	AT7003-100	Label 100 µg Antibody
VcMMAF	AT7003-mIgG1-100	Label 100 µg mouse IgG1 Antibody
VcMMAF	AT7003-500	Label 500 µg Antibody
VcMMAF	AT7003-mIgG1-500	Label 500 µg mouse IgG1 Antibody
VcMMAF	AT7003-1000	Label 1000 µg Antibody
DM4	AT7004-100	Label 100 µg Antibody
DM4	AT7004-mIgG1-100	Label 100 µg mouse IgG1 Antibody
DM4	AT7004-500	Label 500 µg Antibody
DM4	AT7004-mIgG1-500	Label 500 µg mouse IgG1 Antibody
DM4	AT7004-1000	Label 1000 µg Antibody