

ANNEXIN V SELECTION GUIDE

Fluorophore	Laser	Ex (nm)	Em (nm)	Brightness
Annexin V-Alexa Fluor® 350	350	343	441	2
Annexin V-iFluor® 350	350	345	450	3
Annexin V-iFluor® 488	488	491	516	9
Annexin V-FITC	488	491	516	7
Annexin V-Alexa Fluor® 488	488	499	520	8
Annexin V-TRITC	532	544	570	1
Annexin V-Cy3	561	555	569	1
Annexin V-iFluor® 555	561	557	570	8
Annexin V-iFluor® 594	561	587	603	6
Annexin V-Alexa Fluor® 594	561	590	618	4
Annexin V-iFluor® 633	640	640	654	10
Annexin V-iFluor® 647	640	656	670	6
Annexin V-Cy5	640	656	670	6
Annexin V-iFluor® 680	640	684	701	3
Annexin V-Cy5.5	640	684	701	3
Annexin V-iFluor® 700	640	690	713	4
Annexin V-Cy7	750	756	779	10
Annexin V-iFluor® 750	750	757	779	5

* Brightness scale is 1 (dimmiest) to 10 (brightest), as a comparison of all included dyes. Brightness is computed as a function of extinction coefficient, quantum yield and correction factor. Experimental brightness additionally dependent upon DOL of reagent used.





Figure 1. Jurkat cells were treated with 1 μ M staurosporine for 4 hours to induce apoptosis. Following treatment, cells were stained with Annexin V-iFluor® 555 conjugate (Cat No. 20072). Nuclei were labeled with Nuclear Green™ DCS1 (Cat No. 17550). Images were acquired on a confocal microscope.

Untreated Control



Treated with 500 nM SS




Figure 2. Fluorescence images of HeLa cells. Cells were treated with (Right) or without (Left) 500 nM staurosporine (SS) at 37 °C for 4 hours. Cells were then incubated with Annexin V-Cy5.5 conjugate (Red) and measured using fluorescence microscope with a Cy5 filter. Viable cells were labeled with Cellbrite™ Orange (Shown as Green).

For Research Use Only. Not for use in diagnostic or therapeutic procedures. © 2024 AAT Bioquest, Inc. All rights reserved. All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of AAT Bioquest.

